

## Datasheet

### Environmental Glass Bodied pH Electrode Jenway Part Number 924 050



This electrode is designed for measurements at low temperatures or in solutions of low ionic strength. It is also suitable for general purpose pH measurements in accurate Q.C. and research work.

Specifications:

Parameter	Specification
Type	Combination
Sealed/Refillable	Refillable
Fill Solution	3M KCl (025 160)
End Type	Bulb
pH Range	0-11
Temperature Range	0-50°C
Reference Type	HgCl
Junction Type	Annular Ceramic
Body Material	Glass
Dimensions	120x12mm
Cable Length	1 metre (nominal)
Termination	BNC
Isopotential	0±20mV
Slope	59.16mV±3mV @ 25°C
Response Time	95% of final reading in 10 seconds
Stability	<1mV(between 1 & 2 minutes)
Impedance	<300MΩ @25°C

This datasheet is meant for your guidance only.